



# appGENE Total RNA Kit Manual

*Fast, Efficient, and Flexible*



## KIT COMPONENTS

| Product Codes   | ARN015                       | ARN016                         |
|---|------------------------------|--------------------------------|
| Pack size   | 50 preps                     | 250 preps                      |
| ALY buffer <sup>ϕ</sup><br>(volume of 100% β-mercaptoethanol to add) <sup>⊕</sup> | 33ml<br>(330μl) <sup>⊕</sup> | 165ml<br>(1.65ml) <sup>⊕</sup> |
| AW1 buffer <sup>ϕ</sup>   | 35ml                         | 175ml                          |
| AW2 buffer<br>(volume of 96-100% ethanol to add)*                                 | 16.5ml<br>(66ml)*            | 82ml<br>(328ml)*               |
| AEB buffer (for RNA elution)  | 5ml                          | 5x 5ml                         |
| EM reagent (anti-foam)  | 1 ml                         | 5 ml                           |
| Nuclease-free water   | 510μl                        | 2.55ml                         |
| DNase I (RNase-free, lyophilised)   | 2 tubes                      | 10 tubes                       |
| 10X DNase I reaction buffer   | 4.5ml                        | 22.5ml                         |
| appGENE RNA spin columns  | 50 columns                   | 5x 50 columns                  |
| Collection tubes (2ml)  | 50 tubes                     | 5x 50 tubes                    |

### BEFORE STARTING:

1. Mix each buffer well, but do not mix the ALY buffer vigorously as it may foam due to its non-ionic detergent content. (To prevent foaming, add EM Reagent to the ALY buffer during the protocol to a final concentration 3.3% in the lysis mixture – see sample preparation section below).
2. Examine the ALY and AW1 buffers. If sedimentation has occurred, incubate the ALY buffer at 50°C, and the AW1 buffer at 37°C, mixing occasionally until the sediment has dissolved. Cool to room temperature.
3. <sup>⊕</sup> - Add 100% β-mercaptoethanol to the ALY buffer, to a final concentration of 1%. The combined ALY buffer and β-mercaptoethanol will remain stable at 2-8°C for four weeks. When isolating in parts, transfer enough of the ALY buffer for one isolation to a separate RNase-free bottle/tube and add β-mercaptoethanol. Marking the bottle after adding β-mercaptoethanol is recommended.
4. <sup>ϕ</sup> - Protect ALY and AW1 buffers from sunlight.
5. \* - Add appropriate amount of 96-100% ethanol to the AW2 buffer; for details, see the instructions on the bottle label and in the table above. Marking the bottle after adding the alcohol is recommended.

6. DNase I is shipped lyophilized. Each vial contains sufficient amount enzyme for 25 preps. Before using for the first time, reconstitute the DNase I lyophilizate in 255µl of nuclease-free water. Incubate 1 minute at room temperature. Mix carefully by inverting the tube several times. DNase I is sensitive to physical denaturation. Therefore, do not vortex DNase I solution. Divide it into aliquots to avoid excessive freeze-thawing. Do not freeze/thaw more than three times. After reconstitution the DNase I should be kept at -20°C and it is stable for 6 months.
7. Prepare a freezing rack for storage of the eluted RNA.
8. Pre-chill a sterile mortar and pestle.
9. Read the entire protocol bearing in mind the Considerations section to read tips about how to avoid issues like poor RNA integrity due to incorrect storage or the use of inappropriately treated consumables, inefficient elution due to lack of washing, and how to amend the elution volume depending on the concentration of RNA required.

## STORAGE

Store the kit at room temperature (18 - 25°C). Avoid evaporation from buffer bottles by ensuring they are tightly closed before storing. Under proper storage the kit will remain stable for at least 12 months.

## DESCRIPTION

The **appGENE Total RNA kit** provides an easy-to-use and rapid method for the isolation of up to 230µg of total RNA from cultured eukaryotic cells and from fresh or frozen tissues. Multiple samples of up to  $1 \times 10^7$  eukaryotic cells or up to 30 mg of tissue can be processed in parallel in fewer than 20 minutes. Purified RNA can be used in many downstream applications such as RT-PCR, RT-qPCR, Northern blotting, nuclease protection assays, and in vitro translation. The product is intended for research use only.

## SAMPLE PREPARATION

### TISSUE

#### 1.1 Homogenisation via Mortar and Pestle

Put  $\leq 30$ mg tissue previously frozen in liquid nitrogen or dry ice into a pre-chilled, sterile mortar. Using a chilled pestle, carefully, but firmly crush the tissue into smaller pieces and then, into a pulp. Transfer the powder thus obtained into a sterile RNase-free 2ml microcentrifuge tube containing 600µl ALY buffer and 20µl EM reagent and go to step 2 of the protocol. After pulping, a thin, sticky layer may be formed, rather than a powder. If this occurs, add 600µl ALY buffer containing 20µl EM reagent to the mortar and reconstitute the tissue by pipetting and then transfer the lysate into a sterile RNase-free 2ml microcentrifuge tube. Retrieve the tissue remains from the pestle as well.

## 1.2 Homogenisation via Mechanical Homogenizer

Place  $\leq 30\text{mg}$  tissue in a sterile RNase-free 2ml microcentrifuge tube, add 100 $\mu\text{l}$  ALY buffer and 3.3 $\mu\text{l}$  EM reagent and carefully homogenize with a sterile homogenizer tip. After homogenization, retrieve the tissue remains from the knife tip by washing it with 500 $\mu\text{l}$  ALY buffer. Combine the fractions thus obtained and transfer the entire volume to a new sterile, RNase-free 2ml microcentrifuge tube. Go to step 1 of the below protocol.

### CULTURED CELLS

Thaw  $10^4$ - $10^7$  cell suspension or adherent cells, fresh or frozen at 37°C. Centrifuge the cells suspended in growth medium or PBS buffer in a 15ml falcon tube or a 1.5-2ml sterile, RNase-free microcentrifuge tube at 400 x g. If a compact cell pellet is not formed, wash the cells twice with 1ml cold PBS buffer. Add 600 $\mu\text{l}$  ALY buffer containing 20 $\mu\text{l}$  EM reagent. Mix thoroughly by vortexing for 30 sec and subsequent pipetting. In some cases when cells tend to form syncytia (myoblasts) or tight connections (epithelial cells) or they are high in number (approx.  $10^7$  cells), it may be difficult to re-suspend them in ALY buffer. In such cases, pipette carefully, using a  $\geq 1000\mu\text{l}$  pipette tip or a sterile syringe. Do not use filter tips. Transfer everything to a new 2ml sterile, RNase-free microcentrifuge tube. Go to step 1 of the below protocol.

### TOTAL RNA EXTRACTION PROTOCOL

1. Vortex either fragmented tissue / cells (with pre-added 600 $\mu\text{l}$  ALY buffer containing 20 $\mu\text{l}$  EM reagent) for 1 minute.
2. If more isolations are carried out ALY buffer may be prepared for all samples adding EM Reagent to the final concentration 3.3%. E.g. For 10 tissue preps, transfer 6ml of ALY Buffer to an RNase and DNase free tube and add 200 $\mu\text{l}$  EM Reagent and use this buffer for lysis.
3. Centrifuge for 2 minutes at 15,000 x g.
4. Transfer the supernatant into an RNase, DNase free 1.5-2ml sterile, RNase-free microcentrifuge tube.
5. Add 600 $\mu\text{l}$  70% ethanol to the transferred supernatant.
6. Mix well by pipetting or vortexing.
7. Transfer 700 $\mu\text{l}$  of the resulting mixture into an appGENE RNA spin column placed in a collection tube.
8. Centrifuge for 15 sec at  $\geq 8000$  x g.

9. Discard the flow-through and re-use the column, together with the collection tube.
10. Transfer the remaining mixture into the same appGENE RNA spin column and centrifuge at  $\geq 8000 \times g$  for 15 sec.
11. Discard the flow-through and place the appGENE RNA spin column in a new collection tube.

**DNA Removal (OPTIONAL):**

- a. Prewash the appGENE RNA spin column with 500 $\mu$ l AW2 buffer and centrifuge for 1 min at 15,000  $\times g$ .
  - b. For each isolation mix 90 $\mu$ l 10x DNase I reaction buffer and 10 $\mu$ l reconstituted DNase I. Mix by inverting the tube.
  - c. Apply 95 $\mu$ l of the above DNase mixture onto the centre of the appGENE RNA spin column.
  - d. Incubate for 15 minutes at room temperature.
  - e. Add 600 $\mu$ l AW1 buffer and centrifuge for 15 sec at  $\geq 8000 \times g$ .
  - f. Discard the flow-through and re-use the collection tube. **Proceed to step 14.**
12. Add 700 $\mu$ l AW1 buffer and centrifuge for 15 sec at  $\geq 8000 \times g$ .
  13. Discard the flow-through and reuse the collection tube.
  14. Add 500 $\mu$ l AW2 buffer and centrifuge for 30 sec at  $\geq 8000 \times g$ .
  15. Discard the flow-through and reuse the collection tube.
  16. Repeat steps 14 & 15.
  17. Centrifuge for 90 sec at 15,000  $\times g$  (read Considerations, Section 1.3 on washing).
  18. Discard the collection tube and flow-through and carefully transfer the appGENE RNA spin column to a sterile, 1.5ml sterile, RNase-free microcentrifuge (read Considerations, Section 1.4 on elution).
  19. Add 50-100 $\mu$ l elution buffer AEB precisely, onto the centre of the appGENE RNA spin column membrane.
  20. Centrifuge at  $\geq 8000 \times g$  for 1 minute.

21. Remove the appGENE RNA spin column and place the tube with the eluted RNA in a freezing rack. The isolated RNA is ready for use in downstream applications or for storage at -80°C.

## CONSIDERATIONS

### 1.1 Tissue Storage Prior to RNA Extraction

Proper sampling and storing of the biological material prior to RNA isolation is crucial to obtaining a high purity RNA sample. After sampling, the material should either be preserved by deep freezing at -80°C or in liquid nitrogen or stored in RNase inactivating buffers at -20°C. Most tissues must be preserved within 30 minutes of sampling. Tissues rich in RNases (pancreas, liver) must be preserved immediately. When isolating from cell cultures, the best results are achieved with fresh material. If storage is unavoidable, discard the supernatant after centrifugation and freeze the cell pellet at -80°C or in liquid nitrogen.

When isolating from more than the recommended quantity of starting material (>30 mg, >10<sup>7</sup> cells), divide the material into several isolations so that each 30mg (or 10<sup>7</sup> cells) of sample material is isolated with a separate buffer and spin column set. If this quantity is exceeded, the homogenizing column may become clogged and/or the isolated RNA may be of low purity.

### 1.2 RNases

RNases are very active enzymes which do not require any cofactors and are resistant to autoclaving at 121°C for 15 minutes. In order to avoid the degrading effect of the enzymes on the RNA, the following recommendations should be followed:

- Use disposable latex, vinyl or nitrile gloves at all times when working with the RNA. Do not touch any items not designed specifically for RNA work.
- If possible, keep the samples at 2-8°C at all stages of the procedure, including centrifugation. Use decontaminated freezing racks instead of ice in order to avoid RNase contamination. Keeping RNA after elution in the freezing racks is mandatory.
- Disposable plastic ware (tips, tubes) should be RNase-free or autoclaved at 134°C for 18-20 minutes.
- Re-useable plastic ware, glass and porcelain should be soaked overnight in 0.1 N NaOH/0.1% DEPC water (or RNase-free water) and then washed with 0.1% DEPC water (or RNase-free water). When applicable, glass and porcelain (mortars) should be parched at 150-140°C for 2-4 hrs and cooled to room temperature.
- Wipe surfaces, pipettes, centrifuge (wipe the rotor separately) and tube racks with 3% hydrogen peroxide or <0.5% sodium hypochlorite (or any other commercially available RNase inactivating fluid). Prior to decontamination, test the decontaminant on a small area of the material for possible undesired reactions.

### 1.3 Washing

The wash buffer contains alcohol, which may interfere with some enzymatic reactions and also decrease the elution efficiency. It is therefore vital to remove the alcohol completely from the appGENE RNA spin column prior to elution.

### 1.4 Elution

The optimal volume of the AEB elution buffer used should be chosen in line with the amount of the sample material and the final RNA concentration expected. The use of 30-50µl AEB is recommended when extracting from 2-10mg of tissue or <math>10^4</math> cells, as is increasing the elution buffer volume to 100µl when isolating from 10-30mg of tissue or  $10^4$ - $10^7$  cells. If a high RNA concentration is desired, the elution volume may be reduced to 50µl. It should be noted that this may reduce the efficiency of the RNA retrieval. It is essential to apply the elution buffer precisely to the centre of the membrane. When more sample material is to be used for isolation (not recommended as the column can then easily become clogged), full RNA retrieval can be obtained by performing a second elution (100µl). For the second elution, repeat steps 19 to 21 of the above protocol, placing the appGENE RNA spin column in a new, sterile 1.5 ml microcentrifuge tube. AEB buffer does not contain EDTA, which may interfere with some enzymatic reactions.

### 1.5 DNA contamination

All the biological material used for RNA isolation also contains DNA. There is no RNA isolation method which guarantees complete DNA removal unless the RNA sample is treated with DNase after isolation. Even slight DNA contamination (several gDNA copies per reaction) may give an additional signal in a quantitative PCR analysis after reverse transcription. The appGENE Total RNA kit provides efficient on-column digestion of the DNA during RNA purification. The DNase I is removed by wash buffer AW1.

### 1.6 Average DNA Yields from Different Sample Types

| Sample Material   | Amount | DNA Yield |
|-------------------|--------|-----------|
| 293 HEK cell line | $10^6$ | 5 µg      |
| HeLa cell line    | $10^6$ | 7 µg      |
| Liver             | 10mg   | 75 µg     |
| Heart             | 5 mg   | 10 µg     |
| Brain             | 5 mg   | 5 µg      |
| Colon             | 10mg   | 16 µg     |

## TROUBLE SHOOTING / TECHNICAL SUPPORT

For troubleshooting please visit [www.appletonwoods.co.uk/appGENE-Total-RNA-kit-troubleshooting.pdf](http://www.appletonwoods.co.uk/appGENE-Total-RNA-kit-troubleshooting.pdf) for a trouble shooting guide on RT.

## ASSOCIATED PRODUCTS

| Product                              | Pack Size     | Product Code |
|--------------------------------------|---------------|--------------|
| appSCRIPT cDNA synthesis kit         | 25 reactions  | ARP601       |
| appGREEN 1-Step Extreme Low ROX Kit  | 200 reactions | ARP742       |
| appGREEN 1-Step Extreme High ROX Kit | 200 reactions | ARP752       |
| appGREEN 1-Step Opti Low ROX Kit     | 200 reactions | ARP792       |
| appGREEN 1-Step Opti High ROX Kit    | 200 reactions | ARP802       |
| appTAQ Polymerase                    | 500 units     | ARP001       |
| appTAQ RedMix (2X)                   | 200 reactions | ARP062       |
| appTAQ Hot Start Polymerase          | 250 units     | ARP021       |
| appHiFi Polymerase                   | 200 units     | ARP041       |
| Molecular biology grade water        | 100mL         | BMW001       |

More pack sizes available at [www.appletonwoods.co.uk](http://www.appletonwoods.co.uk)