**Instruction Manual**

Corning® LSE™ Digital Dry Bath

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Digitaler, trockener Heizblock Corning LSE

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### Corning® LSE™ Digital Dry Bath

#### Table Of Contents

1.0	Specifications .....	2
2.0	Safety Precautions .....	2
3.0	Installation .....	3
4.0	Controls .....	3
5.0	Operation.....	4
6.0	Calibration.....	4
7.0	RS232 Interface .....	5
8.0	Troubleshooting Guide .....	5
9.0	Technical Support/Service .....	6
10.0	Cleaning and Maintenance .....	6
11.0	Warranty Statement .....	7
12.0	Product Disposal .....	8
13.0	Appendix .....	8

#### Intended Use

Corning LSE single and dual Digital Dry Bath Incubators provide comprehensive designs for a wide variety of life science research applications. Excellent temperature control can deliver accurate and reliable experimental results from one experiment to another.

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

## 1.0 Specifications

Temperature range	Ambient +5°C to 150°C	
Temperature display resolution	0.1°C 4 digit LED	
Temperature uniformity	± 0.2°C (at 37°C in block)	
Temperature accuracy	± 0.3°C	
Temperature controller	Microprocessor – user calibratable	
Timer	1 to 999 minutes in 1 minute increments or continuous	
Operating altitude	Not to exceed 2,000 meters	
Control	PI microprocessor controller	
I/O	RS232 unidirectional	
Dimensions (W x D x H)	7.9 x 10.4 x 3.3 in. (20 x 26.5 x 8.3 cm)	
Weight	4.8 lbs./2.2 kg	
Electrical requirements		
D1100	115V	50/60Hz, 0.97A, fuse 1.6 AT
D1100	230V	50/60Hz, 1.48A, fuse 0.63 AT
D1200	115V	50/60Hz, 1.45A, fuse 2.0 AT
D1200	230V	50/60Hz, 0.72A, fuse 1.0 AT

## 2.0 Safety Precautions

- Do not** use this product in an explosive environment.
- Do not** use in the presence of flammable or combustible material.
- Do not** heat substances that react violently when heated.
- Do not** touch block when hot or when unit is heating (use block lifter).
- Do not** touch area around block or block well when unit is hot.
- Do not** spill liquids into the well area or into the unit side vent holes.
- Connect** unit only to a properly grounded outlet.

**Caution: HOT Surfaces that can burn or cause injury.**

Use of this product in any manner not specified by the manufacturer or modification of the product may cause injury and/or may void the warranty.

## 3.0 Installation

Upon unpacking the Corning® LSE™ Digital Dry Bath, inspect for damages. Shipping damage is the responsibility of the carrier. Inspect that the following are present: instruction manual, power cord, block lifter.

Select a location that is dry and not subject to drafts or moving air from heating or air conditioning vents, or air blown by other equipment. Place the unit on a flat, non-flammable surface. Allow sufficient room around the unit for access and cooling. Six inches minimum on all sides is suggested. Plug the unit into a properly grounded outlet. By using the block lifter, insert the sample block into the well. The unit is now ready for use.

## 4.0 Controls

**On/Off** rocker switch, behind the unit: turns primary power on and off

**Start/Stop** LED, red: starts or stops the unit, illuminates when unit is in heating mode and off in temperature set mode

**Set Temp** LED, red: illuminates when unit is in temperature set mode and off when in heating mode


**Heating** LED, red: illuminates when unit is in heating mode and unit is actually applying heat to the block, on continuously during heat-up and cycles on and off when the unit is at the set temperature

**Timer** LED, green: sets or selects timer mode, illuminates while running in timer mode

**Up** arrow key: raises set temperature or time when unit is in set mode

**Down** arrow key: lowers set temperature or time when unit is in set mode

**Mode** key: switches the unit between heating mode, set temperature mode, or timer mode

## 5.0 Operation

1. Place the Corning® LSE™ Digital Dry Bath on a sturdy and level surface in a safe, dry place away from laboratory traffic.
2. Ensure that the AC power switch is **Off**, and then plug the three-pronged power cord into a grounded three-prong AC outlet of the appropriate voltage (120V or 230V as indicated on the rating sticker near the AC cord on the back of the unit).
3. Select suitable sample block (purchased separately) or appropriate water volume and put it/them into the Corning LSE Digital Dry Bath.
4. Turn the AC power **On**.
5. Run temperature calibration procedure when using the instrument for the first time.
6. Press the **Up** or **Down** arrow keys to adjust temperature.
7. To set the timer press the **Mode** key, and then press the **Up** or **Down** arrow keys to adjust the time.
8. The unit will stop with an alarm when time is up.
9. Press the **Start/Stop** key to start heating.
10. To reset timer during heating, press the **Start/Stop** key to deactivate heating.
11. Press the **Start/Stop** key again to stop the unit.

## 6.0 Calibration

Calibration allows the unit temperature display to be adjusted or matched to the temperature of a single sample or to a calibrated thermometer, making an independent temperature measure of the block.

The Corning LSE Digital Dry Bath is calibrated at the factory at 40°C using a standard small hole-pattern block. If you are using a block with a high heat loss rate, such as a block with large holes or a platform style block, you may choose to recalibrate the unit to your specific application. Also, if you are using very loose or odd shaped vessels, the calibration function can help you match the display temperature to your actual sample temperature.

To calibrate the unit for a given block or sample, first turn the unit off using the **On/Off** switch. Place a thermometer of known accuracy\* into the block thermometer hole, a thermocouple or other sensor into the sample. Make sure there is a good fit between the thermometer and the block or good contact between the sensor and sample or sensor and block.

The unit will then begin to automatically adjust its operating temperature to your original set point with the recalibration factor included. Allow sufficient time for the unit to re-equilibrate and then again compare the thermometer or sensor reading to the calibrated display. They should closely match. If not, repeat the calibration procedure.

## 7.0 RS232 Interface

The Corning® LSE™ Digital Dry Baths have an RS232 unidirectional data port. An optional software disk and interface cable are available that allow a user to use a desktop or laptop computer to record and/or print a record of the temperature profile produced by the dry bath.

## 8.0 Troubleshooting Guide

Symptom	Solution
Display/LEDs do not light up	<ol style="list-style-type: none"><li>1. Check power cord and outlet.</li><li>2. Check <b>On/Off</b> switch.</li><li>3. Check fuse.</li><li>4. Contact Corning for service.</li></ol>
Unit not heating	<ol style="list-style-type: none"><li>1. Check if set point below room temperature.</li><li>2. Check if <b>Start</b> LED illuminated.</li><li>3. Press <b>Start</b> key.</li><li>4. Contact Corning for service.</li></ol>
Unit display overshoots	Normal operation. Display set point during heat-up overshoots on initial heat-up, but block and sample do not overshoot. See Operation section of this manual.
Block or sample temperature not same as display temperature	<ol style="list-style-type: none"><li>1. Check that unit is in heating mode.</li><li>2. Check if unit is sitting in drafty location.</li><li>3. Check accuracy of thermometer.</li><li>4. Check that thermometer is making good contact with the sample block.</li><li>5. Follow calibration procedure.</li></ol>

\*Thermometers used for calibration purposes should have a written calibration certificate and be traceable back to NIST or some other certified body. General lab thermometers are often not accurate enough for calibration work.

## **9.0 Technical Support/Service**

If you have a question about the Corning® LSE™ Digital Dry Baths or have a service inquiry, contact Corning Customer Service immediately at: 1.800.492.1110 (United States and Canada), 1.978.442.2200 (outside the United States), or Contact your local Corning sales office.

Before returning any unit for service, a Return Authorization (RA) number must first be obtained. Equipment sent in without our prior authorization will be returned at the customer's expense. When returning a unit to Corning for service, it should be sent in the original packaging. If this is not possible, be sure that the unit is sufficiently packed. Any damage resulting from improper packaging is the responsibility of the customer. A written explanation should accompany the unit along with the RA number.

## **10.0 Cleaning and Maintenance**

**Make sure that the dry bath and block are cool and the power cord is disconnected before performing any cleaning or maintenance. Repair or maintenance should only be performed by an authorized service technician.**

The dry bath may be cleaned with a moist cloth containing a mild soap solution. Do not immerse the dry bath in water or any liquid.

The blocks may also be cleaned in a mild soapy solution. Be sure that all items have thoroughly dried before attempting to connect the cord or use the unit.

**Spills:** In the event liquid is accidentally spilled into the bath or well area, disconnect the power cord from the outlet and turn the unit upside down to minimize liquid contact with the internal components. Remove the bottom cover and inspect to ensure liquid has not contacted heater elements, electronic controls, or connectors. Have an authorized service technician clean the unit and replace any damaged parts.

## **11.0 Warranty Statement**

Corning Incorporated (Corning) warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. CORNING DISCLAIMS ALL OTHER WARRANTIES WHETHER EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE. Corning's sole obligation shall be to repair or replace, at its option, any product or part thereof that proves defective in material or workmanship within the warranty period, provided the purchaser notifies Corning of any such defect. Corning is not liable for any incidental or consequential damages, commercial loss or any other damages from the use of this product.

This warranty is valid only if the product is used for its intended purpose and within the guidelines specified in the supplied instruction manual. This warranty does not cover damage caused by accident, neglect, misuse, improper service, natural forces, or other causes not arising from defects in original material or workmanship. This warranty does not cover motor brushes, fuses, light bulbs, batteries, or damage to paint or finish. Claims for transit damage should be filed with the transportation carrier.

In the event this product fails within the specified period of time because of a defect in material or workmanship, contact Corning Customer Service at the following numbers: USA: 1.800.492.1110; Canada: 1.978.442.2200. For other regions of the world, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or see the included instruction manual for a list of worldwide support offices.

Corning's Customer Service team will help arrange local service where available, coordinate a return authorization number, or provide a return authorization number and shipping instructions. Products received without proper authorization will be returned. All items returned for service should be sent postage prepaid in the original packaging or other suitable carton, padded to avoid damage. Corning will not be responsible for damage incurred by improper packaging. Corning may elect for onsite service for larger equipment.

Some states do not allow limitation on the length of implied warranties or the exclusion or limitation of incidental or consequential damages. This warranty gives you specific legal rights. You may have other rights which vary from state to state.

No individual may accept for, or on behalf of Corning, any other obligation of liability, or extend the period of this warranty.

For your reference, make a note of the model and serial number, date of purchase, and supplier here.

Model No. \_\_\_\_\_ Serial No. \_\_\_\_\_

Date Purchased \_\_\_\_\_

Supplier \_\_\_\_\_

## 12.0 Product Disposal

According to Directive 2012/19/EU of the European Parliament and Council of 4th July 2012 on waste and electronic equipment (WEEE) as amended, Corning® LSE™ Digital Dry Bath is marked with the crossed-out wheeled bin and must not be disposed of with domestic waste.

Consequently, the buyer shall follow the instructions for reuse and recycling of waste electronic and electrical equipment (WEEE) provided with the products and available at the following link: [www.corning.com/weee](http://www.corning.com/weee).

## Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

The electrical warning indicates the presence of a potential hazard which could result in electrical shock

**ATTENTION:** Hot Surface!

Indicates disposal instruction.

DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. To ensure utmost protection of the global environment and minimize pollution, please recycle this unit.

Register your product warranty online at  
[www.corning.com/lifesciences/warranty](http://www.corning.com/lifesciences/warranty).

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

**Corning Incorporated**  
*Life Sciences*  
836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476  
[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

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For a listing of trademarks, visit [www.corning.com/clstrademarks](http://www.corning.com/clstrademarks).
All other trademarks are the property of their respective owners.
© 2016 Corning Incorporated. All rights reserved. Printed in Mexico 7/16 CLS-EQ-114