TECHNE

N° ICE

OPERATOR'S MANUAL

CONTENTS


	PAGE
INTRODUCTION	2
BEFORE USE	3
Unpacking	5
SAFETY AND INSTALLATION	4
SÉCURITÉ ET CONSIGNES D'INSTALLATION	6
SICHERHEITS- UND INSTALLATIONSINFORMATIONEN	8
INFORMAZIONI SULLA SICUREZZA E L'INSTALLAZIONE	10
INFORMACIÓN DE SEGURIDAD E INSTALACIÓN	12
CONTACT INFORMATION	14
SPECIFICATION	15
OPERATION	16
Installation	16
Power up Screen	16
Setting the Operating Temperature	17
Changing the Temperature	18
After Use	18
User Maintenance	18
FAQS	19
ADDITIONAL INFORMATION	19
Fault Finding	19
Fuses	20
Insulation Tests	20
REPLACEMENT PARTS	20

INTRODUCTION

The N° ICE is a Peltier controlled cooling device ideal for incubating any number or size of sample tubes at sub-ambient temperatures for as long as required. It is an excellent alternative to wet ice buckets which can be messy, inconvenient and lead to sample loss or contamination.

The removable bucket that sits inside the unit is filled with ceramic-coated beads that both maintain the user-set temperature and support the sample vessels. The ceramic beads are chemically resistant and if required, can be autoclaved at temperatures up to 134°C.

The temperature range of the N° ICE (0 to 40°C) means that it is also suitable for other applications such as ligation and enzyme reactions where a constant temperature close to ambient temperatures is required. The unit is supplied with a dark plastic lid which maintains the cooled environment and helps to protect light-sensitive reagents.

BEFORE USE

Before using the N° ICE please make sure you have read this manual carefully. If there is any doubt relating to the proper use of this equipment, the staff at Bibby Scientific Ltd. or your supplier will be happy to assist you.

UNPACKING

When unpacking the unit please ensure that the following have been removed from the packaging:

- N° ICE
- Aluminium bucket
- Bag of ceramic-coated beads
- Mains cables (UK, EU and US)
- Operator's manual
- Guarantee card

The user is advised to keep the original packaging in case the instrument ever needs to be returned for service or repair. Bibby Scientific Ltd. accepts no responsibility for damage incurred unless the unit is correctly packed and transported in its original packaging.

SAFETY AND INSTALLATION

Please read all the information in this manual before using the N° ICE.

WARNING

Contact with low-temperature objects may cause cold burns. Users should be aware of the following potential hazards:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS
- DO NOT place any liquid directly into the instrument.

OPERATOR SAFETY

All operators of Techne equipment must have available the relevant literature needed to ensure their safety. It is important that only suitably trained personnel operate this equipment, in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Bibby Scientific Ltd. the protection provided by the equipment to the operator may be impaired.

All Techne instruments are designed to conform to international safety requirements and are fitted with an over-temperature cut-out. If a safety problem should be encountered, switch off the unit at the mains socket and remove the plug from the electricity supply.

INSTALLATION

The instrument should be carried using both hands. Never move or carry the instrument when in use or connected to the mains electricity supply.

1. All Techne instruments are supplied with a power cable; this may be integral or plug-in.
2.  Before connecting the instrument to the mains electricity supply, check the voltage against the rating plate (located on the back of the unit). Please note that the unit must be earthed to ensure proper electrical safety. Connect the mains cable to a suitable plug according to the table below.

Connection	220/240V, 50/60Hz Supply	110V/120V Supply, 50/60Hz
Live	Brown	Black
Neutral	Blue	White
Earth	Green/yellow	Green

3. UK ONLY: The fused plug supplied with the mains cable is fitted with a 5 amp fuse to protect the instrument and the operator.
4. The units are rated to operate at 100-230V, 50/60Hz.
5. Place the unit on a suitable flat bench or in a fume cupboard if required, ensuring that the air vents on the rear and underside are free from obstruction.

6. Plug the mains cable into the socket on the back of the instrument.
7. Switch on the instrument. Symbols on or near the power switch of the unit have the following meanings:

- I Mains Switch On
- O Mains Switch Off

REPLACEMENT CABLE

Should the mains lead need replacement, a cable of 1mm² of harmonized code H05VV-F connected to an IEC320 plug should be used. **IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN.**

ENVIRONMENTAL CONDITIONS

The N° ICE is designed to operate under the following conditions:

- Indoor use
- Ambient temperature range +5°C to +40°C
- Altitude to 2000m
- Relative humidity not exceeding 80%
- Mains supply fluctuations not exceeding 10%
- Over voltage category II IEC 60364-4-443
- Pollution degree 2

Note: The control specifications are quoted at an ambient temperature of 20°C. The specification may deteriorate outside an ambient temperature of between 10°C and 30°C.

The instrument has been tested for radio frequency interference and is certified under EN61326.

GUARANTEE

The instrument is guaranteed against any defect in material or workmanship for the period specified on the enclosed guarantee card. This period is effective from the date of purchase; within this period all defective parts will be replaced free of charge provided that the defect is not the result of an accident, misuse or negligence. Servicing under this guarantee should be obtained from the supplier of the instrument.

This manual has been prepared for the convenience of Techne's customers and nothing in this manual shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the unit or components.

Notwithstanding the description and specification(s) of the instruments contained in the operator's manual, Techne reserves the right to make such changes as it sees fit to the instruments or to any of the components.

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CONTACT INFORMATION

For technical, sales or servicing information, contact your local Techne dealer or:

UK

Bibby Scientific Limited
Beacon Road, Stone
Staffordshire ST15 0SA
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 810405
e-mail: techne@bibby-scientific.com
www.techne.com

France

Bibby Scientific France SAS
ZI du Rocher Vert - BP 79
77793 Nemours Cedex
France
Tel: +33 1 64 45 13 13
Fax: +33 1 64 45 13 00
e-mail: bsf@bibby-scientific.fr
www.bibby-scientific.com

Middle East

Bibby Scientific Middle East Ltd.
PO Box 27887, Engomi 2433
Nicosia
Cyprus
Tel: +357 22 660 423
Fax: +357 22 660 424
e-mail: sales@bibby-scientificme.com

North and South America

Bibby Scientific US Inc. t/a Techne Inc.
3 Terri Lane, Suite 10
Burlington, N.J. 08016
USA
Toll Free (in NA): 800-225-9243
Tel: +1 609 589 2560
Fax: +1 609 589 2571
e-mail: labproducts@techneusa.com
www.techneusa.com

Italy

Bibby Scientific Italia Srl
Via Alcide de Gasperi 56
20070 Riozzo di Cerro al Lambro (MI)
Italia
Tel: +39 (0)2 98230679
Fax: +39 (0)2 98230211
e-mail: marketing@bibby-scientific.it
www.bibby-scientific.it

Asia

Bibby Scientific Asia
Room 607, Yen Sheng Centre
64 Hoi Yuen Road, Kwun Tong,
Kowloon, Hong Kong
Tel: +852 3583 1581
Fax: +852 3583 1580
e-mail: bibby@bibby-scientificasia.com
www.bibby-scientificasia.com

SPECIFICATION

Temperature

Temperature range	0.0 to 40.0°C
Temperature accuracy	+/-1°C
Temperature units	C or °F
Cooling technology	Peltier
Temperature display resolution	0.1°C
Temperature display	Orange LED, 5 digits

Dimensions

Height	225mm
Width	240mm
Length	190mm
Weight*	4.5kg

*Shipping weight includes bucket and ceramic beads

Power

Power supply	100 to 230V, 50/60Hz
--------------	----------------------

OPERATION


INSTALLATION

See also the sections on **Safety and Installation**.

1. Place the unit on a suitable flat bench or in a fume cupboard if required, ensuring that the air vents on the rear and underside are free from obstruction.
2. Plug the mains cable into the socket on the back of the instrument.
3. Fill the aluminium bucket with the ceramic-coated beads and place the bucket inside the unit.
4. Connect to the mains electricity supply with the plug provided or one wired correctly for the supply. Switch the power ON from the switch located next to the mains input lead. The front display will then light up.

POWER UP SCREEN

When first switched on, the display will show the name of the product and then show the current temperature of the unit. The unit will then proceed to heat or cool to the set temperature.

FRONT PANEL CONTROLS


The front panel consists of two buttons for setting the unit temperature and a 5 digit LED display.

The UP ARROW Button

When this button is pressed, the set temperature is increased.

The DOWN ARROW Button

When this button is pressed, the set temperature is decreased.

SETTING THE OPERATING TEMPERATURE


1. To view the current set temperature, press either the UP or DOWN ARROW button. The set temperature will be displayed and is indicated by a small point in the top left hand corner of the display.
2. To **increase** the set temperature, press and hold the UP ARROW button. The set temperature will momentarily be displayed and then the temperature display will begin to increase. Continue to press the UP ARROW button until the new set temperature is reached then release the button. To fine tune the temperature settings, a single press of the UP or DOWN ARROW keys will increase or decrease the set temperature by 0.1°C.
3. To **decrease** the set temperature, press and hold the DOWN ARROW button. The set temperature will momentarily be displayed and then the temperature display will begin to decrease. Continue to press the DOWN ARROW button until the new set temperature is reached then release the button. To fine tune the temperature settings, a single press of the UP or DOWN ARROW buttons will increase or decrease the set temperature by 0.1°C.
4. Once the set temperature has been changed, the unit will heat or cool to the new set temperature.

Note: The ceramic beads will take some time to equilibrate to the set temperature; therefore it is recommended that the unit is set up in advance to allow time for the required temperature to be reached. To speed up cooling times, the bucket containing the ceramic beads may be stored in a refrigerator prior to use.

Set Point Indicator

New Set Point Indicator


Changing the set temperature will cause the Set Point Indicator to move position while this is being carried out.

CHANGING THE TEMPERATURE UNITS

The temperature of the display can be shown in either °C or °F.

1. To change the temperature units, first turn the power to the unit OFF.
2. Switch the power to the unit ON whilst holding down the UP ARROW button.
3. The temperature unit selection screen will then be displayed. Release the UP ARROW button and then press it again to toggle between °C and °F.

4. Once the required temperature unit has been selected, after a short time the unit will automatically show the current temperature in the selected units.

AFTER USE

After cooling samples, the bucket and beads may be very cold. Allow to warm to ambient conditions before lifting from the unit or wear protective gloves. Avoid handling cold objects with wet hands.

USER MAINTENANCE

Before cleaning the unit, disconnect it from the power supply. The outer case of the N° ICE may be wiped clean with a damp cloth; ethanol, methanol or formaldehyde can also be used. No part of the case or cover should be immersed in the solvents. Do not use aggressive solvents such as acetone or abrasive cleaners. Allow the unit to dry before using.

Before using any cleaning or decontamination method except those recommended here, the responsible person should check with Techne that the proposed method will not damage the equipment.

FAQ'S

Q1 How can I cool the beads more quickly?

A1 **If you intend to use the N° ICE to keep samples cool, then the bucket and beads can be stored in a refrigerator and then inserted into the unit when required.**

Q2 Is it possible to sterilise the beads?

A2 **Yes. Place the beads in a suitable container and autoclave as you would for sterilising other equipment. The beads can be autoclaved at up to 134°C.**

If you require further technical or application assistance please contact Techne at:

E-mail: technehelp@bibby-scientific.com

Phone: +44 (0)1785 810433

Fax: +44 (0)1785 810405

For servicing information please contact:

Bibby Scientific Ltd.

Service Department

Beacon Road

Stone

Staffordshire

ST15 0SA

E-mail: service@bibby-scientific.com

Phone: +44 (0)1785 810475

Fax: +44 (0)1785 810405

We are continually striving to improve our products and software. If you have any comments and suggestions on how we can do things better please send them to us at: techne@bibby-scientific.com

ADDITIONAL INFORMATION

Note that **this equipment should only be dismantled by properly trained personnel**. Removing the outer cover exposes potentially lethal mains voltages. There are no user serviceable parts within this equipment.

FAULT FINDING

Should you have any problems with your N° ICE which cannot be easily remedied, you should contact your supplier and return the unit if necessary. Please include details of the fault observed and remember to return the unit in its original packing. Bibby Scientific Ltd. accepts no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier, giving the full serial number of the unit.

FUSES

If the display on the front panel is not lit, one of the two fuses may have blown. Check that there is no external cause, such as a faulty plug or lead. Check both fuses and replace the faulty fuse with a new one of the correct value (fuse values are given on the label next to the power inlet). Note that **fuses should only be replaced by a qualified electrician.**

The holder for the two fuses is built into the mains input socket. First remove the power cable and then gently prise the fuse drawer open with a flat-bladed screwdriver or similar tool. Each fuse can be removed by using the screwdriver as a lever.

Exchange the faulty fuse in the fuse holder for a working fuse of the correct value. Finally, replace the fuse drawer in the fuse compartment and push the drawer shut. Fuses which blow repeatedly indicate a serious fault and you should return the unit to your supplier for repair.

INSULATION TESTING

This equipment is fitted with RFI suppression circuitry. Any check of the electrical insulation by means of high voltage dielectric testing (for example as in BS EN 61010-1) must be carried out using only a DC voltage.

This unit contains semiconductor components which may be damaged by electric field effects.

REPLACEMENT PARTS

The following replacement parts can be obtained from Bibby Scientific Ltd. or your Techne dealer:

Product Code	Description
FICEBUCKET	N° ICE spare bucket
FTCOOLBE	N° ICE replacement ceramic beads
HH179(S)	UK 230V mains lead with plug
HH180(S)	European mains lead with plug
HFCABLEUS	US mains lead with plug
700738	Fuse Type T5A, 250V (5A)

Bibby Scientific Ltd

Beacon Road, Stone
Staffordshire, ST15 0SA
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 813748
email: sales@bibby-scientific.com
www.bibby-scientific.com

Bibby Scientific France SAS

ZI du Rocher Vert - BP 79
77793 Nemours Cedex
France
Tel: +33 1 64 45 13 13
Fax: +33 1 64 45 13 00
e-mail: bsf@bibby-scientific.fr
www.bibby-scientific.com

Techne Inc.

3 Terri Lane, Suite 10
Burlington, N.J. 08016, USA
Toll Free: 800-225-9243
Tel: +1 609 589 2560
Fax: +1 609 589 2571
e-mail: labproducts@techneusa.com
www.techneusa.com

Bibby Scientific Italia Srl

Via Alcide de Gasperi 56,
20077 Riozzo di Cerro al Lambro,
Milano, Italia
Tel: +39 (0)2 98230679
Fax: +39 (0)2 98230211
e-mail: marketing@bibby-scientific.it
www.bibby-scientific.it

Bibby Scientific Asia

Room 607, Yen Sheng Centre,
64 Hoi Yuen Road, Kwun Tong,
Kowloon, Hong Kong
Tel: +852 3583 1581
Fax: +852 3583 1580
e-mail: bibby@bibby-scientificasia.com
www.bibby-scientificasia.com