

ecoSHIELD™ ECO NITRILE PF 250

1239-45-8 Ethidium bromide 5%	Breakthrough time: 480 minutes	Class 6
1310-73-2 Sodium Hydroxide, 50%	Breakthrough time: 480 minutes	Class 6
50-00-0 Formaldehyde, 35%	Breakthrough time: 480 minutes	Class 6
7664-93-9 Sulphuric acid, 50%	Breakthrough time: 480 minutes	Class 6
79-06-1 Acrylamide 40%	Breakthrough time: 480 minutes	Class 6
7647-01-0 Hydrochloric Acid, 37%	Breakthrough time: 65 minutes	Class 3
67-63-0 (70%) Isopropanol 70%	Breakthrough time: 43 minutes	Class 2
67-63-0 Isopropanol 100%	Breakthrough time: 33 minutes	Class 2
64-17-5 Ethanol, 70%	Breakthrough time: 26 minutes	Class 1
1336-21-6 Ammonium Hydroxide, 25%	Breakthrough time: 15 minutes	Class 1
64-17-5 Ethanol 100%	Breakthrough time: 9 minutes	Class 0
64-19-7 Acetic Acid 100%	Breakthrough time: 4 minutes	Class 0
1330-20-7 Xylene 98,5%	Breakthrough time: 2 minutes	Class 0
109-89-7 Diethylamine 99,5%	Breakthrough time: 0 minutes	Class 0
141-78-6 Ethyl Acetate 100%	Breakthrough time: 0 minutes	Class 0
67-56-1 Methanol 99,9%	Breakthrough time: 0 minutes	Class 0
67-64-1 Acetone 100%	Breakthrough time: 0 minutes	Class 0
67-66-3 Chloroform 99,8%	Breakthrough time: 0 minutes	Class 0
75-05-8 Acetonitrile 99,9%	Breakthrough time: 0 minutes	Class 0