

- ✓ **Stretch Vinyl**
- ✓ **Powder free**
- ✓ **Green**
- ✓ **Smooth**
- ✓ **S-L**
- ✓ **Ambidextrous**

Made of high quality stretch PVC polymer resin coated with Aloe Vera, GN67 provides Aloe Vera's natural benefits such as moisturising hands whilst working and providing relief to cracked or chapped skin.

Attributes

- ✓ **Manufactured from** synthetic emulsion PVC polymer resin with Aloe Vera coating
- ✓ **Latex free** reducing the risk of hypersensitivity and allergic reactions
- ✓ **Soft smooth formulation** for improved feel and fit reducing fatigue during extended wear
- ✓ **Flexible material** for improved feel and fit reducing fatigue during extended wear
- ✓ **Aloe vera coating** will react to body heat moisturising and hydrating skin to help keep hands soft and healthy
- ✓ **Beaded cuff** for added strength and easier donning
- ✓ **Powder free** to reduce potential dust contamination

Regulatory Compliance

- ✓ Complies with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - EN455 parts 1, 2, 3 and 4: Medical Gloves for Single use
 - CE Category: Cat I
- ✓ Tested to:
 - AQL 1.5
 - ASTM F16771-7
- ✓ Visually inspected for holes and minor defects
- ✓ Random samples are filled with 1L of water and examined for leaks
- ✓ Typical force at break: 3.6 Newtons
- ✓ Shelf life of 3 years



Dimensions	S	M	L
Palm Width (mm)	80 ± 10	95 ± 10	110 ± 10
Length (mm)	Min 240	Min 240	Min 240
Thickness (mm) palm	0.09 ± 0.02	0.09 ± 0.02	0.09 ± 0.02
Thickness (mm) finger	0.14 ± 0.02	0.14 ± 0.02	0.14 ± 0.02



Packing


Inner: 100 gloves per dispenser
 Outer: 10 dispensers per case
 Pallet: 84 cases per pallet

Packaging

Outer



Inner




Dimensions

Inner: 230mm L x 120mm W x 63mm H
Outer: 330mm L x 248mm W x 240mm H

Size Colour Guide



Barcodes	Size	Inner Barcode	Outer Barcode
 5 024951 006008	S	5024951006008	5024951006039
	M	5024951006015	5024951006046
	L	5024951006022	5024951006053

Symbols




Lot number



Compliant with the Medical Devices Directive



Latex free



Temperature limitation



Date of manufacture



Non-sterile



Powder free



Recyclable packaging



Use by



Do not re-use



Acceptable quality level



Keep dry



Manufacturer



Keep away from direct sunlight and fluorescent lighting

GN67/JUNE2014



0044 (0) 208 335 3636



0044 (0) 208 335 0300



sales@hpchealthline.co.uk



HPC Healthline Limited, Colwood House, 158 Garth Road, Morden, Surrey, SM4 4LZ, England



www.hpchealthline.co.uk



www.facebook.com/hpchealthline



#hpchealthline