

AURELIA® Robust 9.0®

Examination Gloves



96895	(XS)	ø 100	10 x 100
96896	(S)	ø 100	10 x 100
96897	(M)	ø 100	10 x 100
96898	(L)	ø 100	10 x 100
96899	(XL)	ø 100	10 x 100



Nitrile



Powder Free



Blue



5.5 +/- 0.2 g



Finger-Textured




5 Years

Premium quality Blue Nitrile examination gloves that offer greater tensile strength and protection than both Latex and Vinyl.

Nitrile examination gloves offer extra comfort and flexibility.

- 5mil Powder-Free Nitrile
- 9 Newton
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Latex Free
- AQL 1.5

Sizes & Dimensions



XS	70 ± 10 mm	min. 240 mm
S	80 ± 10 mm	min. 240 mm
M	95 ± 10 mm	min. 240 mm
L	110 ± 10 mm	min. 240 mm
XL	≥ 110 mm	min. 240 mm

Single wall thickness



- ① min. 0.13 mm
- ② min. 0.10 mm
- ③ min. 0.08 mm

Physical properties

	Unaged	Aged
Tensile Strength	min. 14.0 MPa	min. 14.0 MPa
Elongation	min. 500 %	min. 400 %
Strength at Break	min. 9N	min. 9N

Quality Signs & Certificates

89/686/EEC, 93/42/EWG

ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)

EN 455 1-4

EN 420 : 2003 + A 1 : 2009

EN 374 2-3, EN 388

EN 1186: Food Contact Approved (SGS)

