

## Pipettor Compatibility

Brand	Model	Volume Range	Channels	0.1-10µl	20-200µl	1250µl
Biohit	mLine	0.1-3µl, 0.5-10 µl	1	x		
Biohit	mLine	2-20µl, 10-100 µl, 20-200 µl	1		x	
Biohit	mLine	100-1000 µl	1			x
Biohit	mLine	0.5-10 µl	8/12	x		
Biohit	mLine	5-100 µl	8/12		x	
Biohit	Proline	0.1-2.5 µl, 0.5-10 µl	1	x		
Biohit	Proline	2-20 µl, 5-50 µl, 10-100 µl, 20-200 µl	1		x	
Biohit	Proline	100-1000 µl	1			x
Biohit	Proline	0.5-10 µl	8/12	x		
Biohit	Proline	5-50 µl	8/12		x	
Brand	Transferpette	0.1-1µl, 0.5-10 µl	1	x		
Brand	Transferpette	2-20 µl micro	1	x		
Brand	Transferpette	2-20 µl, 5-50 µl, 10-100 µl, 20-200 µl	1		x	
Brand	Transferpette	100-1000 µl	1			x
Brand	Transferpette	0.5-10 µl	8/12	x		
Brand	Transferpette	2-20 µl	8/12	x		
Brand	Transferpette	2.5-25 µl, 5-50 µl, 10-100 µl, 20-200 µl	8/12		x	
Brand	Transferpette	30-300 µl	8/12			
Brand	Transferpette S	0.1-1 µl, 0.1-2.5 µl, 0.5-10 µl	1	x		
Brand	Transferpette S	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Brand	Transferpette S	100-1000 µl	1			x
Brand	Transferpette S	0.5-10 µl	8/12	x		
Brand	Transferpette S	5-50 µl, 10-100 µl, 20-200 µl	8/12		x	
Eppendorf	Research Plus	0.1-2.5 µl, 0.5-10 µl	1	x		
Eppendorf	Research Plus	2-20 µl micro	1	x		
Eppendorf	Research Plus	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Eppendorf	Research Plus	30-300 µl	1			
Eppendorf	Research Plus	100-1000 µl	1			x
Eppendorf	Research Plus	0.5-10 µl	8/12	x		
Eppendorf	Research Plus	10-100 µl	8/12		x	
Eppendorf	Reference	0.1-2.5 µl, 0.5-10 µl	1	x		
Eppendorf	Reference	2-20 µl micro	1	x		
Eppendorf	Reference	2-20 µl, 10-100 µl, 50-200 µl	1		x	
Eppendorf	Reference	100-1000 µl	1			x
Gilson	Pipetman Classic, G, L: P2	0.2-2 µl, 1-10 µl	1	x		
Gilson	Pipetman Classic, G, L: P20	2-20 µl, 20-100 µl, 50-200 µl	1		x	
Gilson	Pipetman Classic, G, L: P1000	200-1000 µl	1			x
Gilson	Pipetman P10L	0.5-10 µl, 2-20 µl, 20-200 µl	8/12		x	

Sales & Customer Service  
**Fax: +44 (0)121 458 5510; Tel: +44 (0)121 458 7740**  
**Email: sales@appletonwoods.co.uk**

www.appleton-scientific.com

## Pipettor Compatibility

Brand	Model	Volume Range	Channels	0.1-10 µl	20-200 µl	1250 µl
HTL	Discovery Comfort & Pro	0.2 - 2.0 µl, 0.5 - 10 µl	1	X		
HTL	Discovery Comfort & Pro	2 - 20 µl	1		X	
HTL	Discovery Comfort	5 - 50 µl	1		X	
HTL	Discovery Comfort & Pro	10 - 100 µl	1		X	
HTL	Discovery Comfort & Pro	20 - 200 µl	1		X	
HTL	Discovery Comfort & Pro	100 - 1000 µl	1			X
HTL	Discovery Comfort	1000 - 10000 µl	1			X
HTL	Discovery Comfort & Pro	1 - 10 µl	8/12	X		
HTL	Discovery Comfort & Pro	5 - 50 µl	8/12		X	
HTL	Discovery Comfort & Pro	20 - 200 µl	8/12		X	
Nichiryo	Nichipet Premium, EX, EX Plus	0.1-2 µl	1	x		
Nichiryo	Nichipet Premium, EX, EX Plus	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Nichiryo	Nichipet Premium, EX, EX Plus	100-1000 µl	1			x
Nichiryo	Nichipet 7000	0.5-10 µl	8/12	x		
Nichiryo	Nichipet 7000	5-50 µl, 40-200 µl	8/12		x	
Rainin	Pipet-Lite SL-2XLS	0.1-2 µl, 0.5-10 µl	1	x		
Rainin	Pipet-Lite SL-20XLS	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Rainin	Pipet-Lite SL-1000XLS	100-1000 µl	1			x
Rainin	Classic PR-2	0.1-2 µl, 0.5-10 µl	1	x		
Rainin	Classic PR-20	2-20 µl, 10-100 µl, 20-200 µl	1		x	
Rainin	Classic PR-1000	100-1000 µl	1			x
Socorex	Calibra® 822	0.2-2 µl, 1-10 µl micro	1	x		
Socorex	Calibra® 822	1-10 µl, 2-20 µl, 10-100 µl, 20-200 µl	1		x	
Socorex	Calibra® 822	100-1000 µl	1			x
Socorex	Calibra® 852	1-10 µl	8	x		
Socorex	Calibra® 852	10-100 µl, 20-200 µl	8/12		x	
Socorex	Acura 825	0.1-2 µl, 0.5-10 µl	1	x		
Socorex	Acura 825	1-10 µl, 2-20 µl, 5-50 µl, 10-100 µl, 20-200 µl	1		x	
Socorex	Acura 825	100-1000 µl	1			x
Socorex	Acura 855	0.5-10 µl	8/12	x		
Socorex	Acura 855	5-50 µl, 20-200 µl	8/12		x	
Thermo/Finn	Finnpette 4500	0.2-2 µl, micro, 0.5-5 µl, micro	1	x		
Thermo/Finn	Finnpette 4500	2-20 µl, 5-50 µl, 10-100 µl, 20-200 µl	1		x	
Thermo/Finn	Finnpette 4500	100-1000 µl	1			x
Thermo/Finn	Finnpette 4500	0.5-10 µl	8/12	x		
Thermo/Finn	Finnpette 4500	5-50 µl	8/12		x	

Sales & Customer Service  
**Fax: +44 (0)121 458 5510; Tel: +44 (0)121 458 7740**  
**Email: sales@appletonwoods.co.uk**

www.appleton-scientific.com