

JB Nova - 5L, 12L, 18L, 26L

NEW

- PERFORMANCE** › Set and Forget™ technology - fast heat-up, reliable temperature control
- INTUITIVE** › Simple to use controls, clear bright display
- QUALITY** › 3 year warranty. UK design and manufacture
- PRACTICAL** › Display lock - disables front panel controls
- CONVENIENT** › Drain tap, tray and lid included
- VALUE** › Advanced dry start and run dry protection - prevents costly service repairs

3
year warranty



Applications

- Education - school / college / university teaching lab
- Healthcare - sample warming
- Industrial - QA/QC testing

Ideal choice for schools and colleges

Specifications



JBN5



3 kg
l: 215 mm
w: 335 mm
h: 200 mm

JBN12



6 kg
l: 380 mm
w: 360 mm
h: 225 mm

JBN18



9 kg
l: 590 mm
w: 335 mm
h: 275 mm

JBN26



9 kg
l: 590 mm
w: 335 mm
h: 275 mm

Capacity	5L	12L	18L	26L
Temperature range	ambient +5 to 95°C			
Temperature display and setting resolution	0.5°C			
Stability (DIN 12876-3 @70 °C)	± 0.5			
Working dimensions (mm)	145 x 290 x 115	315 x 290 x 115	495 x 290 x 115	495 x 290 x 165
User defined calibration points	1	1	1	1
Drain tap	-	●	●	●
Heater power 230V/120V (kW)	0.35 / 0.35	0.8/0.8	1.4/1.05	1.4/1.05
Supply voltage (V)	120 or 230 (50-60Hz)			

Essential accessories (full range available on website)

Flat lid with ring sets

Variable hole diameter to accommodate tall vessels whilst reducing evaporation



Polypropylene spheres

Useful alternative to a lid, minimises evaporation and heat loss whilst allowing easy access to vessels in the bath; particularly useful for tall vessels



Test tube racks

Stainless steel - range of sizes available for 0.5 to 1.5 ml micro-tubes or 10 to 30 mm tubes



Raised shelves

Stainless steel - shelf covers half area of bath. Reversible, allows two shelf depths



Alternative water bath ranges:



SUB Aqua Pro - advanced range

- Temperature range +5 to 99°C
- Stability ± 0.2°C
- 8 bath sizes: 2L, 2L shallow, 5L, 12L, 18L, 26L, 5L & 12L dual bath
- 3 temperature presets
- User defined high temperature alarm and cut out

JB Academy - basic

- Simple and practical, great value
- 3 bath sizes: 5L, 12L, 18L



Call us:

+ 44 (0) 1763 260 811

Email:

sales@grantinstruments.com