

Rotor Compatibility Chart for Corning® Centrifuge Bottles

CORNING

Corning Centrifuge Bottle Compatibility

Corning Cat. No.	Type of Rotor	Rotor Number			
		Beckman Coulter®	Sorvall®	IEC	Fiberlite
431840 431842	6 x 250 mL Fixed Angle	JA-14, JLA-16.250	GSA, SLA-1500	877, 872	15-700-517, 15-700-518
431840 431842	4 x 250 mL Swinging Bucket	JS-7.5, JS-4.3	HS-4, SH-3000	963, 947	n/a
431841 431843	8 x 250 mL Swinging Bucket	n/a	H-2000B	n/a	n/a
431841 431843	6 x 250 mL Swinging Bucket	n/a	n/a	259	n/a
431841 431843	4 x 250 mL Swinging Bucket	GH-3.7, TH-4	H-6000B, RTH-750	219, 218, 224, 268, 245, 243	M4, T4
431841 431843	8 x 250 mL Fixed Angle	n/a	n/a	850S	n/a
431844 431846 431840 431842	6 x 500 mL Fixed angle	JA-10, JLA-10.500	GS-3, SLA-1500	878	15-700-515, 15-700-500, F14S
431845 431847 431841 431843	6 x 1000 mL Swinging Bucket	JS-4.2, JS-3.0	H-6000A	966, 981, 949	C60, P60
431845 431847 431841 431843	4 x 1000 mL Swinging Bucket	JS-4.0, JS-5.2	H-4000	256, 276	P40



Corning Polypropylene Centrifuge Bottles



Corning Polycarbonate Centrifuge Bottles