



Temperature Display with High/Low Alarm RLAA5002 & RLAA5002/1



Factory Fitted
RLAA5002

The **Labcold** battery powered digital Display and Alarm can be factory fitted or purchased as a stand-alone device.

Featuring user adjustable high and low alarm and high and low temperature recording on a clear digital display, it can be factory fitted to the door or panel of the unit depending on cabinet construction.

Made from hygienic plastic, this **Labcold** digital temperature display and alarm is a cost effective solution for daily temperature monitoring and minimum/maximum temperature recording.

Part Number	RLAA5002	RLAA5002/1
Configuration	Factory Fitted	Stand Alone Device
Dimensions (H x W x D)	115 x 62 x 15mm	115 x 62 x 15mm
Temperature Range	-50°C to +70°C	-50°C to +70°C
Resolution	1°C	1°C
Battery	1 AAA (supplied)	1 AAA (supplied)
Digital Temperature Display	✓	✓
Audible Alarm	✓	✓
Alarm Mute	✓	✓
Max/Min Temperature Display	✓	✓
Low Battery Indicator	✓	✓
Detachable Probe	✓	✓



Stand Alone Device
RLAA5002/1