

Haier



APPLETON WOODS LTD



Key Features

- ▶ Most Energy Efficient
- ▶ Low Noise
- ▶ Reliable Performance

More Detail: <http://www.haier.com/uk/commercial/bio-medical/preservation/dw-86l728j.shtml>

Detailed Specifications	
Technical Data	
Type	-86°C ULT Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	HC
Noise Level (dB)	50
Performance	
Cooling performance(°C)	-86
Temperature Range(°C)	-40~-86
Control	
Controller	Microprocessor
Display	LED
Electrical Data	
Power Supply(V/Hz)	220-240/50
Power(W)	1000
Electrical Current(A)	10
Dimensions	
Capacity(L/Cu.Ft)	728/25.7
Net/Gross Weight(approx)	345/385(kg) 760.6/848.8(lbs)
Interior Dimensions(W*D*H)	766*716*1310(mm) 30.2*28.2*51.6(in)
Exterior Dimensions(W*D*H)	1041*980*1980(mm) 36.7*34.5*78.0(in)
Packing Dimensions(W*D*H)	1090*1050*2150(mm) 42.9*41.3*84.6(in)
Container load (20'/40'/40'H)	10/20/20
Functions	
High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High/Low Voltage	N/A
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	3/4
Temperature Recorder	Optional
USB Interface	Y
Rs232/485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional
Racks And Boxes	Optional
CO2 Backup System	
LN2 Backup System	
Chart Recorder	
Others	
Certificate	CE