



# Hotplate Stirrers



# Designed by You



You Select

Model Type – Interface – Surface – Colour

**stuart**<sup>®</sup>

# Built by Us



Ceramic or ceramic-coated aluminium

Temperature sensor compatibility\*

Hot warning light (mains-independent)

Exclusive BioCote® antibacterial treatment

Chemically toughened glass interface

## We Deliver

Your Bespoke Hotplate Stirrer

*\*See specifications for more details*

# Distinctly Different



White\*



Black\*



Blue\*



Pink



Avocado



Tangerine



Red




Ultraviolet

*\*Analogue and Digital models only available in: White, Black & Blue.*

# Analogue Range

With LED display for temperature and stirring



Hotplate Stirrer analogue interface

# Digital Range

With constant digital temperature display



Hotplate Stirrer digital interface

# Premium Range

With TFT display for temperature and stirring



Hotplate Stirrer premium interface

# Technical Specifications

Parameter	Analogue	Digital	Premium
Model types	Hotplate, Stirrer, Hotplate Stirrer	Hotplate, Hotplate Stirrer	
Plate material	Glass Ceramic or Ceramic-coated Aluminium		
Plate dimensions, mm	150 x 150		
Display	Analogue LED temperature and speed arrays	Digital 7 segment temperature display & LED array speed display	Full-colour TFT temperature and speed display
Hot-warning	Mains-independent hot warning light		
Connection for ext. temp sensor	SCT1 contact thermometer*	SCT1 contact thermometer*	Pt100 temperature sensor
Programmable	No	No	Yes
Colour options	White, Black, Blue	White, Black, Blue	White, Black, Blue, Pink, Avocado, Tangerine, Red, Ultraviolet
Heated area, mm	150 x 150 Aluminium, 120 x 120 Ceramic		
Heater power, Watt	700 Aluminium, 500 Ceramic		
Max plate temp, °C	325 Aluminium, 450 Ceramic		
Stirrer speed, rpm	50 - 2000		
Maximum stirring capacity, L	15		
Dimensions (w x d x h), mm	180 x 300 x 93		
Electrical supply	115-230V, 50/60Hz Stirrers 115V or 230V, 50/60Hz Hotplates and Hotplate Stirrers		

\*sold separately



## About Cole-Parmer

Cole-Parmer is a leading global source of Fluid Handling, Test & Measurement and Biosciences products, instrumentation, equipment and supplies. We are proven experts in the fields of temperature measurement, spectroscopy, fluid handling, sample preparation & handling and chromatography. Our portfolio of unique industry-leading brands includes Argos Technologies®, Digi-Sense®, Electrothermal®, Environmental Express®, Ismatec®, Jenway®, Kinesis®, Masterflex®, Oakton®, PCRmax®, Stuart®, Techne®, TELOS®, Traceable®, Vaplock™ and partner with a number of key suppliers with industry leading products.

## Contact Us

### UK

T: +44 (0) 1480 272279  
E: info@coleparmer.co.uk

9 Orion Court, Ambuscade Road,  
Colmworth Business Park, St Neots,  
Cambridgeshire, PE19 8YX

### Germany

T: +49 (0) 9377 92030  
E: sales@coleparmer.de

### France

T: +33 (0) 1486 37800  
E: cpcontact@coleparmer.com

## ColeParmer.com

### Italy

T: +39 (0) 1313 89513  
E: itsales@coleparmer.com

### USA

T: +1 847 549 7600  
E: sales@coleparmer.com