

# Corning® Lambda® Plus Pipettors

CORNING



Lambda Plus single-channel pipettor



Lambda Plus 8-channel and 12-channel pipettors

## Corning Lambda Plus Pipettors

Corning Lambda Plus pipettors have been engineered to provide the highest levels of accuracy, reproducibility, and comfort. Available in single-, eight- and twelve-channel configurations, all Lambda Plus pipettors feature smooth plunger movement and extremely low pipetting forces. Volume can be set by turning either the thumb wheel or the plunger button, and a locking mechanism prevents accidental change.

The Lambda Plus pipettors have been designed to maintain their calibration, even after repeated use and autoclaving. A key is included with each pipettor for easy in-house calibration.

- ▶ Factory calibrated to strict ISO 8655 standards
- ▶ 3-year warranty
- ▶ Fully autoclavable and UV resistant
- ▶ Designed for universal fit with all common brands of pipet tips
- ▶ Adjustable tip ejector

## Corning Lambda Plus Starter Kit

The Lambda Plus starter kit contains four Lambda Plus single-channel pipettors conveniently packaged with useful accessories (manual, calibration key, color-coded rings, tube of piston grease, and identification stickers). The kit includes:

- ▶ Lambda Plus single-channel pipettor, 0.5 to 10  $\mu$ L (Cat. No. 4071)
- ▶ Lambda Plus single-channel pipettor, 2 to 20  $\mu$ L (Cat. No. 4072)
- ▶ Lambda Plus single-channel pipettor, 20 to 200  $\mu$ L (Cat. No. 4074)
- ▶ Lambda Plus single-channel pipettor, 100 to 1000  $\mu$ L (Cat. No. 4075)
- ▶ Universal carousel rack for six single-channel pipettors (Cat. No. 4078)
- ▶ Corning DeckWorks™ pipet tips: Cat. No. 4135 (10  $\mu$ L), Cat. No. 4138 (200  $\mu$ L), and Cat. No. 4140 (1000  $\mu$ L)



Lambda Plus starter kit

## Ordering Information

### Corning® Lambda® Plus Single-channel, 8-channel, and 12-channel Pipettors

Cat. No.	Volume Range (µL)	Channels	Accuracy (%)	Precision (%)	Tip (µL)	Qty/Cs
4070	0.2 - 2.0	1	±12 to ±1.5	<6.0 to <0.7	10	1
4071	0.5 - 10	1	±4.0 to ±0.5	<4.0 to <0.4	10	1
4072	2 - 20	1	±4.0 to ±0.8	<3.0 to <0.4	200	1
4073	10 - 100	1	±1.6 to ±0.8	<0.8 to <0.3	200	1
4074	20 - 200	1	±1.2 to ±0.6	<0.6 to <0.3	200	1
4075	100 - 1,000	1	±1.0 to ±0.6	<0.45 to <0.25	1,000	1
4080	1 - 10	8	±8.0 to ±2.0	<6.0 to <1.2	10	1
4081	5 - 50	8	±4.0 to ±1.6	<2.5 to <0.6	200	1
4082	20 - 200	8	±3.0 to ±1.0	<1.5 to <0.6	200	1
4083	50 - 300	8	±1.6 to ±1.0	<1.5 to <0.6	200/300	1
4084	1 - 10	12	±8.0 to ±2.0	<6.0 to <1.2	10	1
4085	5 - 50	12	±4.0 to ±1.6	<2.5 to <0.6	200	1
4086	20 - 200	12	±3.0 to ±1.0	<1.5 to <0.6	200	1
4087	50 - 300	12	±1.6 to ±1.0	<1.5 to <0.6	200/300	1

### Accessories

4078	Universal carousel rack for six single-channel pipettors
4079	Universal linear rack for six pipettors (six single-channel, maximum of four 8-channel or 12-channel pipettors)
4088	8-channel or 12-channel pipettor stand. Holds one Lambda Plus 8-channel or 12-channel pipettor

### Corning Lambda Plus Starter Kit

Cat. No.	Description	Qty/Cs
4069	Lambda Plus starter kit	1 kit

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, and a complete listing of our International Offices and Distributors, please visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, please call +1.978.442.2200 or contact your local Corning sales office listed below.

#### Corning Incorporated Life Sciences

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476  
[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

#### Worldwide Support Offices

**ASIA/PACIFIC**  
**Australia/New Zealand**  
t 0402-794-347  
**China**  
t 86 21 2215 2888  
f 86 21 6215 2988  
**India**  
t 91 124 4604000  
f 91 124 4604099

**Japan**  
t 81 3-3586 1996  
f 81 3-3586 1291  
**Korea**  
t 82 2-796-9500  
f 82 2-796-9300  
**Singapore**  
t 65 6733-6511  
f 65 6861-2913  
**Taiwan**  
t 886 2-2716-0338  
f 886 2-2516-7500

#### EUROPE

**France**  
t 0800 916 882  
f 0800 918 636  
**Germany**  
t 0800 101 1153  
f 0800 101 2427  
**The Netherlands**  
t 31 20 655 79 28  
f 31 20 659 76 73  
**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

**All Other European  
Countries**  
t 31 (0) 20 659 60 51  
f 31 (0) 20 659 76 73

**LATIN AMERICA**  
**Brasil**  
t (55-11) 3089-7419  
f (55-11) 3167-0700  
**Mexico**  
t (52-81) 8158-8400  
f (52-81) 8313-8589

**CORNING** | **FALCON** **AXYGEN** **GOSELIN** **PYREX**

For a listing of trademarks, visit us at [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks).  
All other trademarks in this document are the property of their respective owners.