

IT CAN'T GET BETTER THAN THIS



iSwix VT (Digital Vortex Mixer)

Speed: 300 to 4200 RPM | Orbital dia. 4 mm



Brushless DC
motor long life



Optional
attachments



Digital display with
speed & time setting



programmable
pulse



3 way switch -
designed for IP 42
class protection



Microprocessor
controlled



Standard Attachment



Round Disc

Optional Attachments



Micro Centrifuge Tube
20 x 1.5/2ml



96 Well Micro Plate



Flask



Cup Attachment

For more details visit us at www.neuation.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.