



APPLETON WOODS LTD

EQUIPMENT & CONSUMABLES
FOR LIFE SCIENCE LABS
www.appletonwoods.co.uk



For the Life Scientist

NEW PCR Plate Centrifuge



Adjustable, compact benchtop PCR plate centrifuge

- Built-in timer and display
- Adjustable speed, up to 2800rpm
- Capacity: 2 x PCR plates, 2 x microtitre plates or 2 x 384 well plates
- Accommodates skirted, non-skirted and semi-skirted plates
- 2 year appleton promise
- Weighs just 6 lbs and has a footprint of 19 x 19.5 cm
- Electronic braking system
- Product code **AES7000**



We are so confident about the quality of this product that it will be **replaced free of charge** should any faults arise within the 2 year warranty period