

## Product Information

### CHO|ONE E

Expression Medium, w/o Insulin, w/o L-Glutamine, with Pluronic™

Cat. No. CHOE-500ML (500 ml), CHOE-1000ML (1000 ml)

### General Information

---

CHO|ONE E allows the cultivation of recombinant CHO cells for fed-batch cultures in spinner flasks or large-scale bioreactors. It can be used for effective cultivation of recombinant CHO cells, achieving high yield of the active protein.

CHO|ONE E already leads to increased cell growth and protein yield but when combined with CHO|ONE Feed 1 and Feed 2 the full potential of the Media System can be deployed.

### Product Specifications

---

Appearance	Clear yellow orange solution
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened store at 4°C and use within 6-8 weeks.
Shipping conditions	Ambient

### Formulation

---

This formulation is our proprietary composition and has no counterparts either in its composition, or in its action.

### Precautions and Disclaimer

---

This product is for research use and further manufacturing only.

### Help Needed?

---

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email ([techservice@capricorn-scientific.com](mailto:techservice@capricorn-scientific.com)) or phone (+49 6424 944640).