

# Getting Started with the Corning® X-SERIES® Platform

CORNING

The Corning X-SERIES Cell Separation Platform helps to simplify your cancer cell therapy workflow by automating many of the necessary steps involved in cell processing.

The Corning X-LAB® and X-WASH® Systems provide a streamlined solution with speed and performance gains isolating mononuclear cells, washing contaminants from cell fractions and concentrating cellular suspensions.

To make it quick and easy to get started, here is a list of all the essential products you'll need to get up and running with the Corning X-LAB and X-WASH Systems.

## Getting Started with the Corning X-LAB System

Automate isolation and harvesting of mononuclear immune cells from blood and blood products with the X-LAB System. Preferred and required products are listed below.

Cat. No.	Description	Qty
6917	Corning X-LAB control module*	1
6918	Corning X-LAB docking station	1
6924	Corning X-LAB disposable cartridges	2
6930	X-SERIES DataTrak software	1
6920	Corning X-Counterweight High	1
6921	Corning X-Counterweight Medium	1
6922	Corning X-Counterweight Low	1
6923	Corning X-Balance ring kit	1

## Getting Started with the Corning X-WASH System

Concentrate cellular suspensions and remove contaminants such as cryoprotectants, spent growth factors, excess viruses, or genetic constructs with the X-WASH System. Preferred and required products are listed below.

Cat. No.	Description	Qty
6928	Corning X-WASH control module*	1
6929	Corning X-WASH docking station	1
6931	Corning X-WASH disposable cartridges	2
6930	X-SERIES DataTrak software	1
6920	Corning X-Counterweight High	1
6921	Corning X-Counterweight Medium	1
6922	Corning X-Counterweight Low	1
6923	Corning X-Balance ring kit	1

\*For cell processing and product functionality convenience a second Control Module is highly recommended. However, only one Control Module is required for operation.



For more specific information on claims, visit [www.corning.com/certificates](http://www.corning.com/certificates).

**Warranty/Disclaimer:** Unless otherwise specified, these products are US class I medical devices and medical device accessories, as general purpose laboratory equipment labeled or promoted for a specific medical use. Corning does not make or support any further claim for these products nor does it conduct any further testing to support any other intended use for the product indicated above. Should you decide that you want to use the product for any intended use other than that made by Corning then you accept that Corning does not make or support that claim and that you as the user are responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of your intended application.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

# CORNING

**Corning Incorporated**  
*Life Sciences*

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

## ASIA/PACIFIC

**Australia/New Zealand**  
t 61 427286832

**Chinese Mainland**  
t 86 21 3338 4338  
f 86 21 3338 4300

**India**  
t 91 124 4604000  
f 91 124 4604099

## Japan

t 81 3-3586 1996  
f 81 3-3586 1291

## Korea

t 82 2-796-9500  
f 82 2-796-9300

## Singapore

t 65 6572-9740  
f 65 6735-2913

## Taiwan

t 886 2-2716-0338  
f 886 2-2516-7500

## EUROPE

CSEurope@corning.com

## France

t 0800 916 882  
f 0800 918 636

## Germany

t 0800 101 1153  
f 0800 101 2427

## The Netherlands

t 020 655 79 28  
f 020 659 76 73

## United Kingdom

t 0800 376 8660  
f 0800 279 1117

## All Other European Countries

t +31 (0) 206 59 60 51  
f +31 (0) 206 59 76 73

## LATIN AMERICA

grupoLA@corning.com

## Brazil

t 55 (11) 3089-7400

## Mexico

t (52-81) 8158-8400