

Freezer Racks

Our freezer racks fit various sizes of freezers and they are compatible with freezer boxes provided by other manufacturers. These freezer racks improve freezer space utilization and satisfy storage requirements of biobanks.

◆Side-access freezer rack features:

- Frame and interior grids are made with 304 stainless steel for durability and better appearance.
- No sharp edges to avoid risk of injury.
- Steel plate thickness: Combination of 0.5 and 1.0 mm.
- Side-access design makes it convenient and economic.
- Suitable for common low temperature freezers.
- Perfectly suitable for standard sized freezer boxes.
- The front side is printed with Haier's stenciled seal.
- Provide personalized customization.



◆Sliding-drawer freezer rack features:

- Made with 304 stainless steel.
- Steel plate thickness: Combination of 0.5 and 1.0 mm.
- Sliding-drawer design for convenient operations.
- Low temperature rated plastic handles are standard on draw racks for ease and safety of operation.
- Labeling slots are standard on draws. The outward pouch design protects the slot from potential damage caused by low temperature.
- The front side is printed with Haier's stenciled seal.
- Provide personalized customization.



Specifications

Product	Model	Rack Description	Rack Dimension (W*D*H)	ULT Freezer Models
2 inch stainless steel rack	DCJ-44-A (2014)	2 inch stainless steel rack 4*4 (2 inch box *16)	138*565*240 (mm) 5.4*22.9*9.4 (in)	DW-338J/490J
	DCJ-44-B (2014)	2 inch stainless steel drawer rack 4*4 (2 inch box *16)	140*565*240 (mm) 5.5*22.2*9.4 (in)	
	DCJ-54-A (2014)	2 inch stainless steel rack 5*4 (2 inch box *20)	138*565*293 (mm) 5.4*22.2*11.5 (in)	DW-388A/486E/ 490/628E
	DCJ-54-B (2014)	2 inch stainless steel drawer rack 5*4 (2 inch box *20)	140*565*293 (mm) 5.5*22.2*11.5 (in)	
	DCJ-55-A (2014)	2 inch stainless steel rack 5*5 (2 inch box *25)	138*685*293 (mm) 5.4*27.0*11.5 (in)	DW-86L388J/418S/418ST/578J/578- JA/578S/578ST/579/579BP/579BPT/728J/728J- A/728S//728ST/729/729BP/729BPT/828J/828JA/ /828S/828ST/829/829BP/829BPT/828W/ DW-86L959BP/959BPT/959W
	DCJ-55-B (2014)	2 inch stainless drawer rack 5*5 (2 inch box *25)	140*685*293 (mm) 5.5*27.0*11.5 (in)	
	DCJ-08-A (2014)	2 inch stainless steel drawer rack 8*1 (2 inch box *8)	140*140*480 (mm) 5.5*5.5*18.9 (in)	
	DCJ-10-A (2014)	2 inch stainless steel rack 10*1 (2 inch box *10)	140*140*600 (mm) 5.5*5.5*23.6 (in)	DW-86W420J/DW-150W200/ DW-150W200J

Haier ULT Freezer + Consumables Solutions

Specifications

Freezer Model	Freezer Rack Model	Description	Freezer Rack Quantity	Freezer Box Quantity	Cryo Tube Quantity (2ml)
DW-86L338J/338JA	DCJ-54-A(2014) (Middle, 2 layers)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	6	120	12000
	DCJ-54-B(2014) (Middle, 2 layers)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	6	120	12000
	DCJ-44-A(2014) (Top and bottom layers)	2 inch stainless steel side-access freezer rack 4*4 (2 inch freezer box*16 units)	6	96	9600
	DCJ-44-B(2014) (Top and bottom layers)	2 inch stainless steel sliding-drawer freezer rack 4*4 (2 inch freezer box*16 units)	6	96	9600
DW-86L388J/418S /418ST/419	DCJ-55- A (2014)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	12	300	30000
	DCJ-55- B (2014)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	12	300	30000
DW-86L486E	DCJ-54- A (2014)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	16	320	32000
	DCJ-54- B (2014)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	16	320	32000
DW-86L490J/490JA	DCJ-54-A(2014) (Top 2 layers)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	8	160	16000
	DCJ-54-B(2014) (Top 2 layers)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	8	160	16000
	DCJ-44-A(2014) (Bottom 2 layers)	2 inch stainless steel side-access freezer rack 4*4 (2 inch freezer box*16 units)	8	128	12800
	DCJ-44-B(2014) (Bottom 2 layers)	2 inch stainless steel sliding-drawer freezer rack 4*4 (2 inch freezer box*16 units)	8	128	12800
DW-86L628E	DCJ-54-A (2014)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	20	400	40000
	DCJ-54-B (2014)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	20	400	40000
DW-86L578J/578JA/ 578S/578ST/579/579BP /579BPT	DCJ-55- A (2014)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	16	400	40000
	DCJ-55- B (2014)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	16	400	40000
DW-86L728J/728JA/ 728S/728ST/729/ 729BP/729BPT	DCJ-55- A (2014)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	20	500	50000
	DCJ-55- B (2014)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	20	500	50000
DW-86L828J/828JA/ 828S/828ST/829/829BP /829BPT/828W	DCJ-55- A (2014)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	24	600	60000
	DCJ-55- B (2014)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	24	600	60000
DW-86L959BP/959BPT /959W	DCJ-55- A (2014)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	28	700	70000
	DCJ-55- B (2014)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	28	700	70000
DW-86L100J	DCJ-53- A (2014)	2 inch stainless steel side-access freezer rack 5*3 (2 inch freezer box*15 units)	4	60	6000
	DCJ-53- B (2014)	2 inch stainless steel sliding-drawer freezer rack 5*3 (2 inch freezer box*15 units)	4	60	6000
DW-86W420J/420JA	DCJ-10- A (2014)	2 inch chest-type freezer rack (2 inch freezer box*10 units)	27	270	27000
DW-86W100J	DCJ-08 - A (2014)	2 inch chest-type freezer rack (2 inch freezer box*8 units)	9	72	72000
DW-150W200	DCJ-10- A (2014)	2 inch chest-type freezer rack (2 inch freezer box*10 units)	12	120	12000
DW-150W200J	DCJ-12-A	2 inch chest-type freezer rack (2 inch freezer box*12 units)	12	144	14400